



Mastercool[®]
"World Class Quality"

**Test Equipment
Depot**
1-800-517-8431

99 Washington Street
Melrose, MA 02176
Phone 781-665-1400
Toll Free 1-800-517-8431

 Visit us at www.TestEquipmentDepot.com

OPERATING INSTRUCTIONS

MULTI-FUNCTION MEDICAL/SURFACE INFRARED THERMOMETER


52225-MED



MULTI-FUNCTION MEDICAL/SURFACE INFRARED THERMOMETER

SPECIFICATIONS:

- Temperature measurement range:
Forehead mode: 93.2 to 108°F (34 to 42.2°C)
Surface mode: -7.6 to 176°F (-22 to 80°C)
- Operating temperature range:
50 to 104°F (10 to 40°C), 15% to 85% RH
- Storage temperature range:
stored at room temperature between -4 to 122 (-20 to 50°C), RH≤85%
Transportation temperature shall be less than 158°F (70°C), RH≤95%
- Atmospheric pressure: 800 to 1013 hPa
- Accuracy:
Forehead mode: ±0.4°F (0.2°C) within 95 to 107.6°F (35 to 42°C)
(Ambient Temp: 59 to 95°F (15 to 35°C), ±0.5°F (0.3°C) for other range
Surface mode: ±0.5°F (0.3°C) within 71.6 to 108°F (22 to 42.2°C), others
±4% or ±4°F (2°C) whichever is greater
- Battery: AAA x 2 pcs
- Battery life: around 3,000 continuous readings
- Expected Service Life: 4 years
- Blue LED Backlight: automatically turned on after measurement, and automatically turned off after 5 seconds
- Enclosure Rating: IP22



NOTE: The thermometer is calibrated at the time of manufacturing. If you question calibration mode, the accuracy of temperature measurements or unexpected events at any time, please contact the service department.

- ⚠ Warning: No modification of this equipment is allowed.
- ⚠ Warning: Do not submerge device into any liquids or expose to direct moisture.




FUNCTIONS

<p>FOREHEAD TEMPERATURE: This thermometer is not meant to replace a visit to the doctor. Please also remember to compare the measurement result to your regular body temperature. Please consult with a doctor if you have health concerns.</p> <ul style="list-style-type: none"> • Please see the “USING THE THERMOMETER” section to learn how to measure the body temperature. 	
<p>SURFACE TEMPERATURE: The surface mode shows the actual and unadjusted surface temperature which is different from the body temperature.</p> <ul style="list-style-type: none"> • Please see the “USING THE THERMOMETER” section to learn how to measure the object temperature. 	
<p>TEMPERATURE INDICATOR: If the thermometer detects a temperature ≥99.5°F (or 37.5°C) under forehead mode, three short beeps will sound followed by one long beep to warn the user that the temperature is above 99.5°F (37.5°C).</p>	
<p>MEMORY LOCATIONS: There are 25 sets of measurement conditions for body temperature.</p> <ul style="list-style-type: none"> • When power is on, press the “ON/MEM” button to see the temperature records with icon. 	

°F / °C: In “power off” mode, press and hold the “START” button, then press the “ON/MEM” button for 3 seconds, the “°F” icon will be switched to the “°C” icon. You can also use the same process to change the LCD display from °C to °F.

NOTE: Memory clear and celsius/fahrenheit switch are together. When you switch celsius and fahrenheit, the memory will be cleared.



MUTE MODE:

The device is set with the sound on, you can set the sound on/off under mute mode.

When power is on, press and hold the “ON/MEM” button for 3 seconds. The 🔊 icon will flash on the LCD screen and then release the “ON/MEM” button to set mute. Thus you will not hear the beep sounds. You can also use the same process to turn off the mute function.

NOTE: If you keep pressing “ON/MEM” button for 2 seconds after 🔊 icon flashing, the device will be powered off WITHOUT setting mute.



USING THE THERMOMETER

NOTE: If there is any temperature difference between the places where the device is stored and where you are going to measure, subject and the device should stay in the same room for at least 15 minutes before measurement.

1. Always make sure the probe lens is clean and without any damage and that the forehead is clean.

⚠ Choking from swallowing small parts and batteries by children or pets is possible, please keep small parts and batteries at places where children and pets cannot reach.

2. Power on:

Press the “ON/MEM” button (see figure 1).



Figure 1

3. Measuring body temperature on the forehead:

Press the “ON/MEM” button to power on the device. **Forehead mode is the default mode.** You can see the 🌡 icon on the screen and you will hear two beep sounds (see figure 1). In this mode, you can hold the thermometer within 1.5” (4 cm) from the central forehead and press the “START” button to get the forehead measurement. The time lapse for measurement might be 1 second. After each forehead measurement, wait for the 🌡 icon to stop flashing to be ready for next measurement.

NOTE:

- Forehead temperature is displayed in oral mode. This mode converts the forehead temperature to display its “oral-equivalent” value.
- Before the measurement, the subject should stay in a stable environment for 5 minutes and avoid exercise/bath for 30 minutes.
- Remember to keep the forehead area clean and away from sweat, cosmetics and scars while taking temperature.
- The “Clinical Bias” is -2.5 ~ -3.1°F (-1.4 ~ -1.7°C).
- The “Limits of Agreement” is 0.98.
- The “Repeatability” is 0.36°F (0.20°C)

MEASURING SURFACE TEMPERATURE

1. After power is on, press and hold the “ON/MEM” button, and press the “START” button one time for “infrared thermometer” mode you will see the 🔍 icon on your LCD display. In this mode, you can get the target surface temperature.

2. When you press the “START” button, you will get the real time temperature immediately. If you press and hold the “START” button, the measurement reading will be continuously updated.

NOTE: This mode shows the actual and unadjusted surface temperature which is different from the body temperature.

POWER OFF

1. The device will automatically shut-off if left idle for more than 1 minute to extend battery life.
2. Manually power off the device by pressing the “ON/MEM” button.

IMPORTANT NOTES

CLEANING AND STORAGE:

Please make sure the probe is clean to ensure an accurate reading.

△ The probe lens is the most delicate part of the thermometer. Use with care when cleaning the probe lens to avoid damage.

- Use alcohol or cotton swabs moistened with 70%~75% alcohol to clean the probe lens.
- Allow the probe to fully dry for at least 1 minute.
- Keep the unit dry and away from any liquids and direct sunlight.
- Storage temperature range: stored at room temperature between -4~122°F (-20~+50°C), RH≤85%
- The probe should not be submerged into liquids.

△ Holding the thermometer too long may cause a higher ambient temperature reading of the probe. This could make the body temperature measurement lower than usual.

BATTERY REPLACEMENT:



When the “low battery” icon indicates the battery is low, the battery should be replaced immediately with AAA *2 pcs batteries.

1. Open the battery cover, use your thumbs to push battery cover out. (see figure 1)
2. Insert the new AAA (*2 pcs) batteries. (see figure 2)
3. Replace the battery cover. (see figure 3)



Figure 1




Figure 2



Figure 3

TROUBLESHOOTING:

Error Message	Problem	Solution
Er	Error 5~9, the system is not functioning properly.	Unload the batteries, wait for 1 minute and re-power. If the message reappears, contact the service department.
Er 1	Measurement before device stabilization.	Wait for “Er1” to disappear.
Er 3	The ambient temperature is not within the range between 50°F~104°F (10°C~40°C).	Allow the thermometer to rest in a room for at least 15 minutes at room temperature: 50°F~104°F (10°C~40°C).
H_i	(1) In forehead mode: temperature taken is higher than 108°F (42.2°C) (2) In surface mode: temperature taken is higher than 176°F (80°C)	Please select the target within specifications. If a malfunction still exists, please contact the service department.
Lo	(1) In forehead mode: temperature taken is lower than 93.2°F (+34°C) (2) In surface mode: temperature taken is lower than -7.6°F (-22°C)	
	Device cannot be powered on to the ready stage.	Change with a new battery.

MANUFACTURER'S DECLARATION-ELECTROMAGNETIC EMISSIONS

Intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user should assure that it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - guidance (for home healthcare environment)
RF emissions CISPR 11	Group 1	This unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.


MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

Intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - guidance (for home healthcare environment)
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical home healthcare environment.

MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

Intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - guidance (for home healthcare environment)
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range. b) interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the thermometer is used exceeds the applicable RF compliance level above, it should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the unit.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE THERMOMETER

Intended for use in an electromagnetic environment (for home healthcare) in which radiated RF disturbances are controlled. The customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER W	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER m		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3 \sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

MANUFACTURER'S DECLARATION-ELECTROMAGNETIC IMMUNITY TEST SPECIFICATIONS FOR ENCLOSURE PORT IMMUNITY TO RF WIRELESS COMMUNICATIONS EQUIPMENT

Intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user should assure that it is used in such an environment.

TEST FREQUENCY (MHz)	BAND a) (MHz)	SERVICE a)	MODULATION b)	MAXIMUM POWER (W)	DISTANCE (m)	IMMUNITY TEST LEVEL (V/m)	COMPLIANCE LEVEL (V/m) (for home healthcare)
385	380 - 390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27	27
450	430 - 470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28	28
710 745 780	704 - 787	LTE Band 13,17	Pulse modulation b) 217 Hz	0,2	0,3	9	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28	28
1720 1845 1970	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28	28
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28	28
5240 5500 5785	5100 - 5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9	9

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50% duty cycle square wave signal.

c) As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

