This manual contains information for how to use and maintain the dissolved oxygen electrode correctly. It also includes explanation for possible problems during operation. Please read it carefully and keep it in a safe place.

Dissolved Oxygen Electrode Instruction Manual

Electrode Model	Order Number	Description
STD011	30031639	Galvanic dissolved oxygen electrode
Maintenance kits	30693682	Maintains kits STD011

Specifications


Electrode type: galvanic dissolved oxygen electrode
 Measuring range: 0-20 mg/L (or ppm)
 Measuring range: 0-200 % sat
 Response time: 60 seconds to reach 95% of the final reading
 Temperature range: 0-50°C
 Cable length: 1 m
 Electrode length: 120 mm
 Electrode diameter: 12 mm
 Cable interface: BNC

Introduction

STD011 dissolved oxygen electrode uses galvanic principle to measure dissolved oxygen concentration in water or solution. It provides a millivolt signal proportional to the oxygen concentration in the water. This electrode cannot measure oxygen directly. The concentration (mg/L or PPM) is calculated according to the dissolved oxygen, temperature, salinity and total atmospheric pressure. The dissolved oxygen electrode is composed of cathode, anode and electrolyte. The electrolyte is separated from the sample through oxygen permeation membrane. The oxygen passing through the film reacts with the cathode and releases electrons to produce an electric current.

To use with a temperature electrode, you can choose STTEMP30 model (83032959).

Electrode Assembly

1. Unscrew the protective bottle.
2. Unscrew the membrane cap from the top of the electrode. Do not touch the membrane because the oil on the skin will damage it. Fill the membrane cap with electrode filling solution.
3. Put the membrane cap carefully back on the electrode. During this process, the filling

solution will overflow.


4. Rinse the electrode and place it in distilled water. **Note:** Do not completely submerge the membrane cap.
5. After the assembly is complete, connect the electrode BNC interface to the meter. Use it after warming up and polarizing for at least 10 minutes for the first time.

Zero Calibration

1. After the electrode is warmed up and polarized, start the calibration procedure according to the meter's instrument manual.
2. Take the electrode out of the water, shake to dry it, and put it in the zero oxygen reagent (sodium sulfite calibration solution).
3. When the display reading is stable, the program will stop automatically.

Air (Water-saturated Air) Calibration

1. Rinse the electrode with distilled water and gently shake it to dry.
2. Add water to the bottle to a depth of about 1 cm and screw the lid and probe back onto the bottle as shown in the graphic.
Important note: during this step, do not touch the film on the front of the electrode.

3. When the display reading is stable, the program will stop automatically.
4. The face icon of the meter indicates the calibration slope:

Slope: 80-125%, which is close to the theoretical value. The electrode is in good condition.

Slope: 60-80%, which deviates from the theoretical value. Please clean the electrode or replace the electrode membrane cap.

Measurement

1. Gently stir the electrode in the water. Monitor the dissolved oxygen concentration in

real-time readings.

2. Please keep stirring the sample. During measurement, water must always flow through the electrode due to oxygen consumption. If the electrode measures a calm sample, the dissolved oxygen reading will continue to decrease, resulting in instability and incorrect readings.
3. If the water is not at room temperature, in order to ensure that the temperature electrode reaches equilibrium, let the electrode stay in the sample for at least two minutes, and turn on the automatic temperature compensation function.
4. After calibration and measurement, unscrew the electrode BNC interface and separate it from the meter. Clean the electrode head with distilled water, and screw the protective bottle onto it.

Short-term Storage (Less Than 24 Hours)

Place the electrode in a moist electrode protection bottle. Keep the membrane cap moist and add distilled water in time.

Long-term Storage (More Than 24 Hours)

Remove the membrane cap and rinse the inside and outside of the cap with distilled water. Shake the membrane cap dry. Also rinse and dry the exposed internal components of the anode and cathode (absorb dry with laboratory wipes). Reinstall the membrane cap loosely on the electrode body for storage. Don't tighten it.

Replace the Membrane Cap

The membrane cap needs to be replaced regularly. Replace it once a year if it is not used frequently, or after about 6 months of continuous use. When your dissolved oxygen electrode responds slowly during calibration or measurement readings, we suggest you to replace the membrane cap.

Troubleshooting

If the electrode calibration data is not ideal or the reading deviation is large, the following points need to be considered:

- Has the protective bottle been removed?
- Is there a filling solution in the electrode cap? If so is the solution free of air bubbles?
- Has the dissolved oxygen electrode been warmed up?
- Has the dissolved oxygen electrode been calibrated?

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