



- Up to 5,000rpm with automatic ref conversion
  - Safety lid-lock Prevents opening the lid until the rotor comes to a complete stop
  - Unique air flow design reduces increase in sample temperature
  - 8 tube capacity Elaht position rotor accepts tubes ranging from 5ml to 15ml







15ml Blood Collection Tubes



10ml Blood Collection Tubes



5-7ml Blood Collection Tubes\*

Benchmark's LC-8 series of laboratory centrifuges is uniquely designed to meet a wide range of applications in both the clinical and research laboratory. An eight position rotor is included and designed to directly accept 10ml and 15ml testtubes (conical or round bottom), while optional adapters are also available for use with all popular sizes of blood collection and culture tubes. Unlike most compact centrifuges in this class, the rotor can be easily removed from the chamber for cleaning.

**Centrifuge Tubes** 

The back-lit, liquid crystal display provides full control over the set run time and rotational speed. The rotational speed can be set/viewed in rpm or rdand acceleration/deceleration levels can be easily adjusted in the advanced settings mode. To ensure safety, the digital microprocessor ensures that the lid has been closed and that the safety lid lock is engaged prior to starting a centrifugation cycle. In addition, an advanced imbalance detection system constantly monitors the balance and vibration, automatically bringing the rotor to a stop if an imbalance is detected.

## Technical Data:

LC-8 LC-8 Plus 400 - 3500rpm (1500xg) 400 - 5000rpm (3070xg) Speed: Capacity: 8x15ml 8×15ml Accel/Decel: 9 sec. / 25 sec. 19 sec./39 sec. Increment: 100 rpm increments 20 sec. to 99 min. Momentary: Yes **Dimensions**: 11.6 x 14.4 x 8.3 in. 29.5 x 36.5 x 21 cm Weight 15 lbs / 6.8 kg Electrical: 115V, 50Hz or 230V, 60Hz

## Ordering Information:

C3100\* LC-8 Centrifuge with 8 x 15ml rotor (maximum speed: 3500rpm / 1500xg)

LC-8 *Plus* Centrifuge with 8 x 15ml rotor C3200

(maximum speed: 5000rpm / 3070xg)

C3100-ADP Optional adapters for 5ml & 7ml blood collection/culture tubes, 8/pk

C2570 Optional adapters for 5ml centrifuge tubes, 4/pk

\*115V with US plug. To order in 230V, add (-E) to the Item number.





99 Washington Street Melrose, MA 02176 Phone 781-665-1400 Toll Free 1-800-517-8431



Visit us at www.TestEquipmentDepot.com