



OIML weight sets include a complete group of test weights that can meet requirements for testing to OIML standards. E1 weight sets offer the highest level of accuracy for standards of reference. E2 weight sets provide a high level of accuracy and are ideal for calibrating analytical and precision balances. F1 weight sets are used to calibrate precision and legal-for-trade class II balances. M1 weight sets offer a lower level of accuracy, but are affordable and ideal for calibrating compact portable balances and bench scales.

Other classes and value of weights are available upon request.

Stainless steel weight sets for calibration that meets testing requirements and standards

| Model | Class | Density | Weights Included | Accessories* |
|--|---------|------------------------|---|--------------|
| E1 1mg - 50g Calibration Weight Set | OIML E1 | 7.9 g/cm ³ | 1mg, 2mg(2), 5mg, 10mg, 20mg(2), 50mg, 100mg, 200mg(2), 500mg, 1g, 2g(2), 5g, 10g, 20g(2), 50g | 1, 2, 3, 4 |
| E1 1mg - 100g Calibration Weight Set | OIML E1 | 7.9 g/cm ³ | 1mg, 2mg(2), 5mg, 10mg, 20mg(2), 50mg, 100mg, 200mg(2), 500mg, 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g | 1, 2, 3, 4 |
| E1 1mg - 200g Calibration Weight Set | OIML E1 | 7.9 g/cm ³ | 1mg, 2mg(2), 5mg, 10mg, 20mg(2), 50mg, 100mg, 200mg(2), 500mg, 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g | 1, 2, 3, 4 |
| E2 1mg - 200g Calibration Weight Set | OIML E2 | 7.96 g/cm ³ | 1mg, 2mg(2), 5mg, 10mg, 20mg(2), 50mg, 100mg, 200mg(2), 500mg, 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g | 1, 2, 3, 4 |
| E2 1g - 200g Calibration Weight Set | OIML E2 | 7.96 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g | 1, 2, 3, 4 |
| F1 1mg - 200g Calibration Weight Set | OIML F1 | 7.9 g/cm ³ | 1mg, 2mg(2), 5mg, 10mg, 20mg(2), 50mg, 100mg, 200mg(2), 500mg, 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g | 3, 4 |
| F1 1mg - 500g Calibration Weight Set | OIML F1 | 7.9 g/cm ³ | 1mg, 2mg(2), 5mg, 10mg, 20mg(2), 50mg, 100mg, 200mg(2), 500mg, 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g(2), 500g | 3, 4 |
| F1 1g - 100g Calibration Weight Set | OIML F1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g | 4 |
| F1 1g - 200g Calibration Weight Set | OIML F1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g | 4 |
| F1 1g - 500g Calibration Weight Set | OIML F1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g(2), 500g | 4 |
| F1 1g - 2kg Calibration Weight Set | OIML F1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g(2), 500g, 1kg, 2kg | 4 |
| M1 1g - 100g Calibration Weight Set | OIML M1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g | 1, 2 |
| M1 1g - 200g Calibration Weight Set | OIML M1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g | 1, 2 |
| M1 1g - 500g Calibration Weight Set | OIML M1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g(2), 500g | 1, 2 |
| M1 1g - 2kg Calibration Weight Set | OIML M1 | 7.9 g/cm ³ | 1g, 2g(2), 5g, 10g, 20g(2), 50g, 100g, 200g(2), 500g, 1kg, 2kg | 1, 2 |
| *Accessories Included: 1. Brush 2. Wipe Cloth 3. Tweezer 4. Cotton Glove | | | | |

Features

- Stainless steel weights for durability and accuracy
- Aluminium case to safely store weights

About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Supplied by