## Western Blot Boxes BLOTTING CONTAINERS



- A Perfect Box for Every Blot
- Reduce Antibody Usage
- Hinged Lid Minimizes Evaporation
- For Halves, Wholes, Strips, or Minis

B1206 Multibox™ 6

85mm X 85mm (for Novex™)

B1300-8BK

No matter what size your Western blot is, there's a box that fits it perfectly. With boxes that are correctly matched to common blot sizes, you can fit more blots on shakers/rockers while keeping reagent and antibody usage to an absolute minimum. Additionally, all comers are formed at a 90° angle and the bottom surface is uniformly flat, further reducing consumption of costly reagents and antibodies. (For example, popular mini-gels require only about 4 ml.)

At a price that makes them truly disposable, all sizes are egineered from crystal clear polystyrene (except B1200-7BK & B1200-15BK, which are an opeque black that protects light sensitive samples). All sizes are also equipped with a easy-to-open hinged lid that reduces evaporation.

| Cat No.   | Description      | Applications   | Volume   |
|-----------|------------------|--|----------|
| B1200-3   | Micro-strip      | Various  | 1 ml     |
| B1200-6   | Mini-rectangle   | Half mini-protein gels                                       | 3 mi     |
| B1200-7   | Medium Rectangle | Mini-protein gels, Bio-Rad                                   | 4 ml     |
| B1200-9   | Long Strip       | Proteen 3 (cut up)   | 2 ml     |
| B1200-11  | Medium Strip     | Verlous  | 3 mi     |
| B1200-13  | Large            | Carribres: PAGIIr Gold, Mighty<br>Small, Novex SureLock Mini | 8 mi     |
| B1200-16  | Square           | Owl Panguin P808/Puffin                                      | 7 ml     |
| B1200-17  | Micha Large      | Bio-Red Offerion, Ameraham                                   | 11 mi    |
| B1200-18  | XXII.            | Bio-Red Offerion, American                                   | 12 ml    |
| B1206     | Long Strip (x6)  | Protes 3 (multiple)  | 2 ml x 6 |
| B1300-88K | Medium Squere    | Novex Mini Geia  | 4 ml     |



## Ordering Information

| Cat. No.  | <u>Dimensions</u>                                  |
|-----------|--|
| B1200-3   | 2 % x 1% x ¾ in<br>(7.3 x 3.0 x 1.9 cm)            |
| B1200-5   | 2 1/2 x 2 x 1 1/4 in<br>(7.3 x 5.1 x 3.2 cm)       |
| B1200-7   | 3 1/2 x 2 1/18 x 1 in<br>(8.9 x 6.5 x 2.5 cm)      |
| B1200-7BK | 3 1/2 x 2 1/16 x 1 in<br>(8.9 x 6.5 x 2.5 cm)      |
| B1300-8BK | 8.6 x 8.8 x 2.8cm<br>(pull off lid)                |
| B1200-9   | 3 3/4 x 1 3/18 x 5/8 ln<br>(9.5 x 3.0 x 1.6 cm)    |
| B1200-11  | 4 3/8 x 2 1/18 x 1 1/8 in<br>(11.1 x 5.2 x 4.4 cm) |

| Cat. No.   | Dimensions                                   |
|------------|--|
| B1200-13   | 4 % x 3 ½ x 1 % in<br>(11.7 x 8.9 x 2.8 cm)  |
| B1200-13BK | 4 % x 3 ½ x 1 % in<br>(11.7 x 8.9 x 2.8 cm)  |
| B1200-15   | 4 ½ x 4 ½ x 1 ¼ in<br>(11.4 x 11.4 x 2.8 cm) |
| B1200-15BK | 4 ½ x 1 ½ x 1 ¼ in<br>(11.4 x 11.4 x 2.8 cm) |
| B1200-17   | 6 x 4 x 1 1/4 ln<br>(15.2 x 10.2 x 2.8 cm)   |
| B1200- 19  | 6 x 4 x 2 in<br>(15.2 x 10.2 x 5.1 cm)       |
| B1206      | 6 compartment box<br>(11.7 x 8.9 x 2.8 cm)   |